

Submittal Data Sheet



Complete assembly of pressure/vacuum switch, 4" diameter pressure/vacuum gauge, and gas specific DISS union demand check valve. Entire assembly is cleaned for oxygen service.

PX-PSM-24 shown above

PSM Assembly Features

Unique to Powerex, the PSM assembly consists of a Tee Assembly using the Powerex Genesys[™] manifold high/low line pressure transducer and DISS union demand check valve – both of which are included with the manifold and eliminate the need to purchase an additional high/low source line pressure switch – with the 4" diameter pressure gauge.

PST / PSX Pressure Switch Specifications

Pressure and vacuum switches for the medical gas alarm system shall incorporate a U.L. listed single pole double throw snap-action micro switch. Normally open or normally closed circuitry shall be field-selected to be compatible with the medical gas alarm system. Electrical rating shall be 10 amperes at 120 volt AC.

PST / PSX Switch Adjustment

Low pressure gases such as oxygen, nitrous oxide and medical air are typically set to alarm at 40 psig (low line pressure) and 60 psig (high line pressure). High pressure nitrogen are typically set to alarm at 160 psig (low line pressure) and 200 psig (high line pressure). The vacuum switch is pre-set at the factory to alarm at 12 in/Hg. Switches are field adjustable (if necessary) by following instructions included.

Ordering Information

| Part # | Gas Service and Pressure Switch Used |
|-------------|---------------------------------------------------------------------------------------------------|
| PX-PST-04 | N2O with High/Low 1-80 psig switch and 0-100 psi gauge |
| PX-PST-08 | CO2 with High/Low 1-80 psig switch and 0-100 psi gauge |
| PX-PSX-12 | N2 with two single setting pressure switches (high and low) and 0-300 psi gauge |
| PX-PST-16 | MED AIR with High/Low 1-80 psig and 0-100 psi gauge |
| PX-PST-22 | VAC with single setting 1-30" Hg switch and 0-30" Hg gauge |
| PX-PST-24 | OXY with High/Low 1-80 psig switch and 0-100 psi gauge |
| PX-PST-HYP | Hyperbaric Oxygen with High (84)/Low (56) switch and 0-100 psi gauge |
| PX-PSX-2080 | Instrument Air with High (200) and Low (160) switches and 0 – 300 psi gauge |
| PX-PST-32 | WAGD with single setting 1-30" Hg switch and 0-30" Hg gauge |
| PX-PSM-04 | N2O assembly with transducer and DISS union provided as part of Genesys™ manifold |
| PX-PSM-08 | CO2 assembly with transducer and DISS union provided as part of Genesys [™] manifold |
| PX-PSM-12 | N2 assembly with transducer and DISS union provided as part of Genesys [™] manifold |
| PX-PSM-16 | MED AIR assembly with transducer and DISS union provided as part of Genesys [™] manifold |
| PX-PSM-24 | OXY assembly with transducer and DISS union provided as part of Genesys [™] manifold |
| PX-PSM-2080 | Instrument Air assembly with transducer and DISS union provided as part of Genesys™ manifold |



PX-PSX-12 shown above

PST / PSX Pressure Switch Features

- Removable copper diaphragm capsule
- U.L. and C.S.A. listed.
- Instrument quality snap-action switch
- Normally open or normally closed contacts



Pressure Switch Wiring in PST / PSX assemblies

Use this chart for United Electric (U. E.) Hi/Low Pressure Switches Only

| Circuit | Common | Normally Closed | Normally Open |
|--------------------|--------|-----------------|---------------|
| Low line pressure | * | | * |
| High line pressure | * | * | |

* Indicates the electrician should connect wires to these terminals.

Use this chart for United Electric (U. E.) Vacuum Switch Only

| Circuit | Common | Normally Closed | Normally Open |
|------------|--------|-----------------|---------------|
| Low vacuum | * | | * |

* Indicates the electrician should connect wires to these terminals.

Use this chart for Barksdale Hi/Low Pressure Switches Only

| Circuit | Common | Normally Closed | Normally Open |
|--------------------|---------|-----------------|---------------|
| Low line pressure | *Purple | | *Red |
| High line pressure | *Brown | *Orange | |

* Indicates the electrician should connect wires to these terminals.

Use this chart for Barksdale Vacuum Switch Only

| Circuit | Common | Normally Closed | Normally Open |
|------------|---------|-----------------|---------------|
| Low vacuum | *Purple | *Red | |

* Indicates the electrician should connect wires to these terminals.