


Adapters, Couplers & Fittings





For more than half a century, Allied Healthcare Products has been supplying quality adapters, couplers and fittings to our customers. Our history dates back to the medical department of National Cylinder Gas (NCG), which produced gases in cylinders for hospitals in the 1930s and early 1940s. In 1953, NCG installed its first piped medical gas system in a hospital. Over the years the NCG Medical Division evolved into the Chemetron Medical Corporation, and finally became Allied Healthcare Products, Inc. in 1980. In 1986, Allied acquired Oxequip Healthcare Industries, a major manufacturer of construction and patient care products since 1936.

Through such acquisitions and transformations we have sought to strengthen our leadership position in medical gas, respiratory therapy, and emergency medical products throughout the world.

In an age when the world has become a global village, we at Allied Healthcare Products continually renew our commitment to the "Global Support of Life."

Chemetron® Style Quick-Connect Adapters

1/8" NPT MALE

Oxygen	12-01-1000
Air (rectangular striker)	12-01-1100
Vacuum (rectangular striker)	12-01-1200
Nitrous Oxide	12-01-1300
WAGD	12-01-1400
CO ₂	12-01-1500
<i>For Gas Outlets 1966 and Earlier:</i>	
Air (round striker)	12-01-1101
Vacuum (round striker)	12-01-1201



1/8" NPT FEMALE

Oxygen	12-01-2000
Air (rectangular striker)	12-01-2100
Vacuum (rectangular striker)	12-01-2200
Nitrous Oxide	12-01-2300
WAGD	12-01-2400
<i>For Gas Outlets 1966 and Earlier:</i>	
Air (round striker)	12-01-2101
Vacuum (round striker)	12-01-2201



▶	Oxygen	12-11-1000	Air	12-11-1100
	Vacuum	12-11-1200	Nitrous Oxide	12-11-1300

1/4" NPT MALE

Oxygen	12-01-1010
Air	12-01-1110
Vacuum*	12-01-1210
Nitrous Oxide	12-01-1310
WAGD	12-01-1410
CO ₂	12-01-1510



1/4" HOSE BARB

Oxygen	12-01-5000
Air (rectangular striker)	12-01-5100
Vacuum (rectangular striker)	12-01-5200
Nitrous Oxide	12-01-5300
WAGD	12-01-5400
CO ₂	12-01-5500
<i>For Gas Outlets 1966 and Earlier:</i>	
Air (round striker)	12-01-5101
Vacuum (round striker)	12-01-5201



DISS MALE

Oxygen	12-01-3000
Air	12-01-3100
Vacuum	12-01-3200
Nitrous Oxide	12-01-3300



5/16" HOSE BARB

Vacuum	12-01-5202
WAGD	12-01-5402



DISS FEMALE HEX NUT

Oxygen	12-01-4001
Air	12-01-4101
Vacuum	12-01-4201
Nitrous Oxide	12-01-4301



DISS FEMALE HAND TIGHT

Oxygen	12-01-4000
Air	12-01-4100
Vacuum	12-01-4200
Nitrous Oxide	12-01-4300
CO ₂	12-01-4500



DISS KNURLED NUT

Air	12-01-4102
Vacuum	12-01-4202
Nitrous Oxide	12-01-4302



Chemetron® Style Quick-Connect Couplers

QUICK-CONNECT TO 1/8" MALE

Oxygen	11-02-0021
Vacuum	11-02-0026
Air	11-02-0031
WAGD	11-02-0035
Nitrous Oxide	11-02-0046
CO ₂	11-02-0047



QUICK-CONNECT TO 1/4" NPT MALE

Oxygen	11-02-0023
Vacuum	11-02-0025
Nitrous Oxide	11-02-0028
Air	11-02-0030
CO ₂	11-02-0050
WAGD	11-02-0052



QUICK-CONNECT TO 1/4" HOSE BARB

Oxygen	11-02-0020
Vacuum*	11-02-0024
Nitrous Oxide	11-02-0027
Air	11-02-0029
CO ₂	11-02-0032
WAGD	11-02-0034



Oxygen	11-02-1023	Air	11-02-1030
Vacuum	11-02-1025		

QUICK-CONNECT TO DISS MALE

Oxygen	11-02-0022
Vacuum	11-02-0041
Nitrous Oxide	11-02-0042
Air	11-02-0043
WAGD (CGA 2220)	11-02-0044
CO ₂	11-02-0045



Oxygen	11-02-1020	Air	11-02-1029
Vacuum	11-02-1024		

QUICK-CONNECT TO 5/16" HOSE BARB

Vacuum	11-02-0037
WAGD	11-02-0039



Oxygen	11-02-5022	Air	11-02-5043
Vacuum	11-02-5041	CO ₂	11-02-5045

Chemetron® Quick-Connect Duplex Couplers

QUICK-CONNECT TO DOUBLE COUPLER CENTER

Oxygen	11-01-0010
Vacuum	11-01-0015
Air	11-01-0027



QUICK-CONNECT TO DOUBLE COUPLER OFFSET LEFT

Oxygen	11-01-0011
Vacuum	11-01-0112
Air	11-01-0113



QUICK-CONNECT TO DOUBLE COUPLER OFFSET RIGHT

Oxygen	11-01-0012
Vacuum	11-01-0115
Air	11-01-0111



QUICK-CONNECT TO DOUBLE COUPLER — "Y" STYLE

Oxygen	11-03-0001
Vacuum	11-03-0003
Air	11-03-0002




*Recommend use of 5/16" hose minimum for vacuum service.

Ohmeda® Style Quick-Connect Adapters

1/8" NPT MALE


Oxygen	12-02-1000
Air	12-02-1100
Vacuum	12-02-1200
Nitrous Oxide	12-02-1300
WAGD	12-02-1400
CO ₂	12-02-1500



1/8" NPT FEMALE

Oxygen	12-02-2000
Air	12-02-2100
Vacuum	12-02-2200
Nitrous Oxide	12-02-2300
WAGD	12-02-2400



	Oxygen	12-12-1000	Air	12-12-1100
	Vacuum	12-12-1200	Nitrous Oxide	12-12-1300

1/4" NPT MALE

Oxygen	12-02-1010
Air	12-02-1110
Vacuum	12-02-1210
Nitrous Oxide	12-02-1310
WAGD	12-02-1410
CO ₂	12-02-1510



1/4" HOSE BARB

Oxygen	12-02-5000
Air	12-02-5100
Vacuum	12-02-5200
Nitrous Oxide	12-02-5300
WAGD	12-02-5400
CO ₂	12-02-5500



DISS MALE

Oxygen	12-02-3000
Air	12-02-3100
Vacuum	12-02-3200
Nitrous Oxide	12-02-3300



DISS FEMALE HAND TIGHT

Oxygen	12-02-4000
Air	12-02-4100
Vacuum	12-02-4200
Nitrous Oxide	12-02-4300
CO ₂	12-02-4500



DISS FEMALE HEX NUT

Oxygen	12-02-4001
Air	12-02-4101
Vacuum	12-02-4201
Nitrous Oxide	12-02-4301
CO ₂	12-02-4501



*New style Evacuation (CGA 2220) after 1994.
For old style (pre-1994), contact Allied's Service/Repair Department.

Ohmeda® Style Quick-Connect Couplers

QUICK-CONNECT TO 1/4" NPT MALE

Oxygen	11-10-0042
Vacuum	11-10-0044
Nitrous Oxide	11-10-0043
Air	11-10-0045



QUICK-CONNECT TO 1/4" HOSE BARB

Oxygen	11-10-0050
Vacuum	11-10-0051
Nitrous Oxide	11-10-0052
Air	11-10-0053



QUICK-CONNECT TO 1/8" NPT MALE

Oxygen	11-10-0039
Vacuum	11-10-0046
Nitrous Oxide	11-10-0048
Air	11-10-0049



QUICK-CONNECT TO 1/8" NPT FEMALE

Oxygen	11-10-0055
Vacuum	11-10-0056
Nitrous Oxide	11-10-0057
Air	11-10-0058



QUICK-CONNECT TO DISS MALE

Oxygen	11-10-0059
Vacuum	11-10-0060
Nitrous Oxide	11-10-0064
Air	11-10-0069



Ohmeda® Style Quick-Connect Duplex Couplers

QUICK-CONNECT TO DOUBLE COUPLER — "Y" STYLE

Oxygen	11-03-0004
Air	11-03-0005
Vacuum	11-03-0006



QUICK-CONNECT TO DOUBLE COUPLER — BAR STYLE

Oxygen	11-10-0090
Vacuum	11-03-0026



Puritan-Bennett Style Quick-Connect Adapters

1/8" NPT MALE

Oxygen	12-03-1000
Air	12-03-1100
Vacuum	12-03-1200
Nitrous Oxide	12-03-1300



1/8" NPT FEMALE

Oxygen	12-03-2000
Air	12-03-2100
Vacuum	12-03-2200
Nitrous Oxide	12-03-2300



1/4" NPT MALE

Oxygen	12-03-1010
Air	12-03-1110
Vacuum	12-03-1210
Nitrous Oxide	12-03-1310



1/4" HOSE BARB

Oxygen	12-03-5000
Air	12-03-5100
Vacuum	12-03-5200
Nitrous Oxide	12-03-5300



DISS MALE

Oxygen	12-03-3000
Air	12-03-3100
Vacuum	12-03-3200
Nitrous Oxide	12-03-3300




DISS FEMALE HAND TIGHT

Oxygen	12-03-4000
Air	12-03-4100
Vacuum	12-03-4200
Nitrous Oxide	12-03-4300



DISS FEMALE HEX NUT

Oxygen	12-03-4001
Air	12-03-4101
Vacuum	12-03-4201
Nitrous Oxide	12-03-4301



▶	Oxygen	12-13-4000	Air	12-13-4100
	Vacuum	12-13-4200	Nitrous Oxide	12-13-4300

Puritan-Bennett Style Quick-Connect Couplers

QUICK-CONNECT TO 1/4" NPT MALE

Oxygen	11-10-0014
Vacuum	11-10-0016
Nitrous Oxide	11-10-0018
Air	11-10-0020



QUICK-CONNECT TO DOUBLE COUPLER — "Y" STYLE

Oxygen	11-03-0013
Vacuum	11-03-0015



Schrader Style Quick-Connect Adapters

1/8" NPT MALE (Non-Swivel)

Oxygen	12-06-1000
Air	12-06-1100
Vacuum	12-06-1200
Nitrous Oxide	12-06-1300




DISS FEMALE HAND TIGHT (Non-Swivel)

Oxygen	12-06-4000
Air	12-06-4100
Vacuum	12-06-4200
Nitrous Oxide	12-06-4300




1/4" NPT MALE (Non-Swivel)

Oxygen	12-06-1010
Air	12-06-1110
Vacuum	12-06-1210
Nitrous Oxide	12-06-1310



1/4" HOSE BARB (Non-Swivel)

Oxygen	12-06-5000
Air	12-06-5100
Vacuum	12-06-5200
Nitrous Oxide	12-06-5300



DISS MALE (Non-Swivel)

Oxygen	12-06-3000
Air	12-06-3100
Vacuum	12-06-3200
Nitrous Oxide	12-06-3300



1/4" HOSE BARB (Swivel)

Oxygen	12-06-5001
Air	12-06-5101
Vacuum	12-06-5201
Nitrous Oxide	12-06-5301



Schrader Style Quick-Connect Couplers

1/4" NPT MALE

Oxygen	11-10-0036
Vacuum	11-10-0033
Nitrous Oxide	11-10-0031
Air	11-10-0035



1/8" NPT FEMALE

Vacuum	11-10-0032
Air	11-10-0034



3/8" NPT FEMALE

Oxygen	11-10-0028
Nitrous Oxide	11-10-0030

Oxequip® Twist Style Quick-Connect Adapters


1/8" NPT MALE

Oxygen	12-05-1000
Air	12-05-1100
Vacuum	12-05-1200
Nitrous Oxide	12-05-1300
CO ₂	12-05-1500
Nitrogen	12-05-1700



1/8" NPT FEMALE

Oxygen	12-05-2000
Air	12-05-2100
Vacuum	12-05-2200
Nitrous Oxide	12-05-2300
CO ₂	12-05-2500
Nitrogen	12-05-2700



1/4" NPT MALE

Oxygen	12-05-1010
Air	12-05-1110
Vacuum	12-05-1210
Nitrous Oxide	12-05-1310
CO ₂	12-05-1510
Nitrogen	12-05-1710



1/4" HOSE BARB

Oxygen	12-05-5000
Air	12-05-5100
Vacuum	12-05-5200
Nitrous Oxide	12-05-5300
CO ₂	12-05-5500



DISS MALE

Oxygen	12-05-3000
Air	12-05-3100
Vacuum	12-05-3200
Nitrous Oxide	12-05-3300
CO ₂	12-05-3500
Nitrogen	12-05-3700
Evacuation (CGA 2220)	12-05-3410



DISS FEMALE HAND TIGHT

Oxygen	12-05-4000
Air	12-05-4100
Vacuum	12-05-4200
Nitrous Oxide	12-05-4300
CO ₂	12-05-4500
Nitrogen	12-05-4700



▶	Oxygen	12-15-4000	Air	12-15-4100
	Vacuum	12-15-4200	Nitrous Oxide	12-15-4300

Oxequip® Medstar® Style Quick-Connect Adapters

1/8" NPT MALE

Oxygen	12-04-1000
Air	12-04-1100
Vacuum	12-04-1200
Nitrous Oxide	12-04-1300
WAGD	12-04-1400
CO ₂	12-04-1500



1/8" NPT FEMALE

Oxygen	12-04-2000
Air	12-04-2100
Vacuum	12-04-2200
Nitrous Oxide	12-04-2300
WAGD	12-04-2400
CO ₂	12-04-2500



Oxygen	12-14-1000	Air	12-14-1100
Vacuum	12-14-1200	Nitrous Oxide	12-14-1300

1/4" NPT MALE

Oxygen	12-04-1010
Air	12-04-1110
Vacuum	12-04-1210
Nitrous Oxide	12-04-1310
WAGD	12-04-1410
CO ₂	12-04-1510



1/4" HOSE BARB

Oxygen	12-04-5000
Air	12-04-5100
Vacuum	12-04-5200
Nitrous Oxide	12-04-5300
WAGD	12-04-5400
CO ₂	12-04-5500



DISS MALE

Oxygen	12-04-3000
Air	12-04-3100
Vacuum	12-04-3200
Nitrous Oxide	12-04-3300
WAGD (CGA 2220)	12-04-3410
CO ₂	12-04-3500



DISS FEMALE HAND TIGHT

Oxygen	12-04-4000
Air	12-04-4100
Vacuum	12-04-4200
Nitrous Oxide	12-04-4300

Oxequip® Medstar® Style Quick-Connect Couplers

DISS MALE

Oxygen	2027-JD
Nitrous Oxide	2027-JC
Air	2027-JA
Vacuum	2027-JE
WAGD (CGA 2220)	2027-HW
CO ₂	2027-JG
HPN	2027-JB



1/4" NPT MALE

Oxygen	2027-HU
Nitrous Oxide	2027-HT
Air	2027-HR
Vacuum	2027-HV



Oxygen	2037-HU	Air	2037-HR
Vacuum	2037-HV	Nitrous Oxide	2037-HT

Oxygen	2037-JD	Air	2037-JA
Vacuum	2037-JE	Nitrous Oxide	2037-JC

QUICK-CONNECT TO DOUBLE COUPLER — "Y" STYLE

Oxygen	11-03-0016
Vacuum	11-03-0018
Air	11-03-0017



Oxygen	11-13-0016	Air	11-13-0017
Vacuum	11-13-0018		



DISS Male Fittings

DISS MALE TO 1/4" HOSE BARB

WITHOUT CHECK

Oxygen	12-80-3005
Air	12-80-3105
Vacuum	12-80-3205
Nitrous Oxide	12-80-3305
Nitrogen	12-80-3705



DISS MALE TO 1/4" HOSE BARB

WITH CHECK

Oxygen	11-30-0006
Nitrous Oxide	11-30-0007
Vacuum	11-30-0008
Air	11-30-0009
CO ₂	11-30-0010
Nitrogen	11-30-0012



DISS MALE TO 1/8" NPT MALE

WITHOUT CHECK

Oxygen	12-80-3001
Air	12-80-3101
Vacuum	12-80-3201
Nitrous Oxide	12-80-3301
WAGD (CGA 2220)	12-80-3418
CO ₂	12-80-3501
Nitrogen	12-80-3701



DISS MALE TO 1/8" NPT MALE

WITH CHECK

Oxygen	12-80-3011
Air	12-80-3111
Vacuum	12-80-3211
Nitrous Oxide	12-80-3311
WAGD (CGA 2220)	12-80-3408
CO ₂	12-80-3511
Nitrogen	12-80-3711



DISS MALE TO 1/4" NPT MALE

WITH CHECK

Oxygen	12-80-3010
Air	12-80-3110
Vacuum	12-80-3210
Nitrous Oxide	12-80-3310
CO ₂	12-80-3510
Nitrogen	12-80-3710



DISS Female Fittings

DISS FEMALE TO 1/4" HOSE BARB

FEMALE HEX NUT

Oxygen	12-80-4005
Air	12-80-4105
Vacuum	12-80-4205
Nitrous Oxide	12-80-4305
CO ₂	12-80-4505
Nitrogen	12-80-4705



Oxygen, green plastic (50/bag) 65401

DISS FEMALE TO 1/4" HOSE BARB

FEMALE KNURLED NUT

Oxygen	12-80-4015
Air	12-80-4115
Vacuum	12-80-4215
Nitrous Oxide	12-80-4315
Nitrogen	12-80-4715



Oxygen, white plastic
(each) 12-80-4045



DISS Female Fittings *(continued)*

DISS FEMALE TO 1/4" HOSE BARB *(continued)*

FEMALE HAND TIGHT NUT

Oxygen	12-80-4035
Air	12-80-4135
Vacuum	12-80-4235
Nitrous Oxide	12-80-4335
CO ₂	12-80-4535
Nitrogen	12-80-4735



DISS FEMALE TO 5/16" HOSE BARB

FEMALE HEX NUT

Vacuum	12-80-4206
WAGD (CGA 2220)	12-80-4429

FEMALE HAND TIGHT NUT

Vacuum	12-80-4207
--------	------------



DISS FEMALE TO 3/16" HOSE BARB

Oxygen Fitting Only	12-90-4000
Oxygen with Chain and Locking Ring	12-90-4001



DISS FEMALE TO 1/8" NPT MALE

HEX NUT

Oxygen	12-80-4021
Air	12-80-4121
Vacuum	12-80-4221
Nitrous Oxide	12-80-4321
WAGD (CGA 2220)	12-80-4427
CO ₂	12-80-4521
Nitrogen	12-80-4721



KNURLED NUT

Oxygen	12-80-4011
Air	12-80-4111
Vacuum	12-80-4211
Nitrous Oxide	12-80-4311
WAGD (CGA 2220)	12-80-4417
CO ₂	12-80-4511
Nitrogen	12-80-4711



DISS FEMALE TO 1/8" NPT MALE *(continued)*

HAND TIGHT NUT

Oxygen	12-80-4031
Air	12-80-4131
Vacuum	12-80-4231
Nitrous Oxide	12-80-4331
WAGD (CGA 2220)	12-80-4437
CO ₂	12-80-4531
Nitrogen	12-80-4731



Oxygen	12-81-4031	Vacuum	12-81-4231
Air	12-81-4131		

DISS FEMALE TO 1/8" NPT FEMALE

KNURLED NUT

Oxygen	12-80-4012
Air	12-80-4112
Vacuum	12-80-4212
Nitrous Oxide	12-80-4312
WAGD (CGA 2220)	12-80-4418
CO ₂	12-80-4512
Nitrogen	12-80-4712




HAND TIGHT NUT

Oxygen	12-80-4032
Air	12-80-4132
Vacuum	12-80-4232
Nitrous Oxide	12-80-4332
CO ₂	12-80-4532
Nitrogen	12-80-4732



DISS FEMALE TO DISS FEMALE ADAPTERS

Oxygen	12-90-4024
Oxygen (90° Elbow)	12-90-4124



DISS MALE TO DISS MALE ADAPTERS

Oxygen	Z205
--------	------




Double Bar & "Y" Block Assemblies

CHEMETRON® QUICK-CONNECT TO DOUBLE DISS MISCELLANEOUS


CHEMETRON® QUICK-CONNECT TO DOUBLE DISS FEMALE CENTER

Oxygen	11-01-0004
Vacuum	11-01-0070
Air	11-01-0049



CHEMETRON® QUICK-CONNECT TO DOUBLE DISS FEMALE OFFSET LEFT

Oxygen	11-01-0005
Vacuum	11-01-0071
Air	11-01-0072



CHEMETRON® QUICK-CONNECT TO DOUBLE DISS FEMALE OFFSET RIGHT

Oxygen	11-01-0006
Vacuum	11-01-0073
Air	11-01-0074



SCHRADER QUICK-CONNECT TO DISS MISCELLANEOUS (Non-SWIVEL)

SCHRADER QUICK-CONNECT TO DISS FEMALE 90° ELBOW

Oxygen	11-10-0007
Vacuum	11-10-0085
Air	11-10-0086



SCHRADER QUICK-CONNECT TO DOUBLE DISS FEMALE CENTER

Oxygen	11-10-0008
Vacuum	11-10-0087
Air	11-10-0088



DISS TO DOUBLE DISS MISCELLANEOUS


ADAPTER DISS MALE TO DOUBLE DISS FEMALE

Oxygen	11-01-0037
Vacuum	11-01-0075
Air	11-01-0076



ADAPTER DISS FEMALE TO DOUBLE DISS FEMALE

Oxygen	11-01-0009
Vacuum	11-01-0109
Air	11-01-0077



ADAPTER DISS FEMALE TO DOUBLE DISS MALE WITH CHECK

Oxygen	11-01-0116
Vacuum	11-01-0117



ADAPTER DISS FEMALE TO DOUBLE DISS FEMALE — "Y" STYLE CONNECTOR

Oxygen	11-03-0007
Vacuum	11-03-0009
Air	11-03-0008



ADAPTER DISS FEMALE TO DOUBLE DISS MALE WITH CHECK — "Y" STYLE CONNECTOR

Oxygen	11-03-0010
Vacuum	11-03-0012
Air	11-03-0011



BLANK DOUBLE ADAPTER BAR

1/8" NPT Female Center x Double	
1/8" NPT Female	11-01-0055



1/8" NPT Female Offset Left/Right x Double 1/8" NPT Female	11-01-0056
--	------------



BLANK "Y" STYLE CONNECTOR


1/8" NPT Female x Double	
1/8" NPT Female	12-90-0028
1/8" NPT Female x Double	
1/4" NPT Female	12-90-0029



NPT Pipe Fittings

NPT CROSS

1/8" Female	12-90-0013
1/4" Female	12-90-0014




NPT DROP EAR

Elbow: 1/8" Female	12-90-0015
Tee: 1/8" Female	12-90-0016



NPT 90° ELBOW

1/8" Female x 1/8" Female	12-90-0001
1/4" Female x 1/4" Female	12-90-0002



NPT 90° STREET ELBOW

1/4" Female x 1/4" Male	12-90-0003
1/8" Female x 1/8" Male	12-90-0004
1/8" Male x 1/8" Male	T12749



NPT 45° STREET ELBOW

1/8" Male x 1/8" Female	12-90-0005
-------------------------	------------



NPT FEMALE TEE

1/8" Chrome	12-90-0006
1/4" Chrome	12-90-0007
1/4" Nickel	12-90-0008



NPT STREET TEE

1/8" Female x Male	12-90-0009
1/4" Female x Male	12-90-0010



NPT PIPE BRANCH TEE

1/8" Female x Male	12-90-0011
1/4" Female x Male	12-90-0012



HEX NIPPLE

1/8" NPT Male x 1/8" NPT Male	12-90-1021
1/4" NPT Male x 1/8" NPT Male	12-90-1121
1/4" NPT Male x 1/4" NPT Male	12-90-1120



NPT TO HOSE BARB

1/8" NPT Male x 3/16" Diameter	12-90-1007
1/8" NPT Male x 1/4" Diameter	12-90-1005
1/8" NPT Male x 5/16" Diameter	12-90-1006
1/4" NPT Male x 1/4" Diameter	12-90-1105
3/8" NPT Male x 1/4" Diameter	12-90-1008
1/8" NPT Female x 1/4" Diameter	12-90-2005
1/4" NPT Female x 1/4" Diameter	12-90-2105
1/4" NPT Female x 5/16" Diameter	12-90-2106



NPT CONNECTORS

1/8" NPT Female x 1/8" NPT Male	12-90-2001
1/8" NPT Female x 1/4" NPT Male	12-90-2000
1/4" NPT Female x 1/8" NPT Male	12-90-2101
1/4" NPT Female 1/4" NPT Male	12-90-2100
1/8" NPT Female x 1/8" NPT Female	12-90-2002
1/4" NPT Female x 1/8" NPT Female	12-90-2102
1/4" NPT Female x 1/4" NPT Female	12-90-2112



PIPE NIPPLE


1/8" NPT Male x 1/8" NPT Male x 1 1/2" Long	12-90-1101
1/8" NPT Male x 1/8" NPT Male x 2" Long	12-90-1201
1/8" NPT Male x 1/8" NPT Male x 3 1/2" Long	12-90-1301



Bottle Holder Adapters

CHEMETRON® QUICK-CONNECT TO SINGLE OUTLET PORT

Vacuum (rectangular striker) 11-01-0016



CHEMETRON® QUICK-CONNECT TO DOUBLE OUTLET PORT

Vacuum (rectangular striker) 11-01-0017



CHEMETRON® DISS FEMALE TO SINGLE OUTLET PORT

Vacuum 11-30-0002



CHEMETRON® DISS FEMALE TO DOUBLE OUTLET PORT

Vacuum 11-30-0001



OHMEDA® QUICK-CONNECT TO SINGLE OUTLET PORT

Vacuum 11-10-0047

Dynacon Nitrogen Fittings

COUPLER, DYNACON TO 1/4" HOSE BARB

Nitrogen 11-20-0012



CHECK UNIT, DYNACON WITH 1/4" NPT MALE FOR CONVERSION OF REGULATORS

Nitrogen 11-20-0008



REPLACEMENT KNOB FOR DYNACON OUTLETS AND COUPLERS

Nitrogen 11-20-0011




ADAPTER, DYNACON MALE TO 1/4" HOSE BARB

Nitrogen 11-20-0001



ADAPTER, DYNACON MALE WITH SET SCREW LOCK ACCEPTING SCHRADER MALE

Nitrogen
Miscellaneous 11-20-0009



Miscellaneous Fittings & Valves

MASTER-FLO

1/8" NPT Female to
1/8" NPT Female

12-90-2052



YOKE BLOCK INDEX TO 1/4" HOSE BARB

Oxygen


11-02-0038



EXTENSION MANIFOLDS (DISS)

	OXYGEN	AIR	VACUUM
2 outlet	12-90-0020	12-90-0030	12-90-0033
3 outlet	12-90-0021	12-90-0031	12-90-0034
4 outlet	12-90-0022	12-90-0032	12-90-0035

Extension manifolds also available with Quick-Connect couplers.



AUXILIARY

Oxygen (9/16"-18" Male
Check)

12-90-0017

Air (3/4"-16" Male Check)

12-90-0018



FERRULE

Many of the adapters, couplers and fittings are used in making hose assemblies.



Although the preparation of hose assemblies by the user is discouraged, information on ferrules is included for convenience.

For all 1/4" hose assemblies (kit of 10) 15-90-0050

For all 5/16" hose assemblies (kit of 10) 15-90-0049

TAMPER-PROOF LOCK-OUT ADAPTERS FOR CHEMETRON QUICK-CONNECT OUTLETS

Oxygen 11-01-0035

Vacuum 11-01-0036

Air 11-01-0041

Additional ISO Styles

BRITISH PROBE TO 1/8" NPT MALE

Oxygen	12-30-1000	Air (MA) (4 Bar)	12-30-1100
Vacuum	12-30-1200	Nitrous Oxide	12-30-1300

GERMAN DIN PROBE TO 1/8" NPT MALE

Oxygen	12-30-1400	Air (MA)	12-30-1500
Vacuum	12-30-1600	Nitrous Oxide	12-30-1700

FRENCH PROBE TO 1/8" NPT MALE

Oxygen	12-30-1800	Air	12-30-1900
Vacuum	12-30-2000	Nitrous Oxide	12-30-2100

Conductive Medical Grade Hose *(UL Listed)*

Chemetron® Bulk Hose is made from carbon impregnated neoprene to provide a conductive path from end fitting to end fitting when hose assembly does not exceed 30 feet in length. Chemetron® Conductive Hose meets all requirements as specified in NFPA 99. This hose is color coded to conform to CGA C-9 standards.

HOSE

	<u>500 Foot Roll</u>	<u>By The Foot</u>
1/4" Green for Oxygen Service	15-09-0501	15-09-0001
1/4" White for Vacuum Service (ISO Oxygen for Canada)	15-09-0502	15-09-0002
5/16" White for Vacuum Service (ISO Oxygen for Canada)	15-09-0506	15-09-0006
1/4" Blue for Nitrous Oxide Service	15-09-0503	15-09-0003
1/4" Yellow for Air Service (ISO Vacuum for Canada)	15-09-0504	15-09-0004
1/4" Black for All Other Services	15-09-0510	15-09-0010

HOSE RETRACTOR

Designed to be installed on ceiling-mounted medical gas outlets but can also be attached to the hose itself. Hose assemblies are kept overhead by the retractor and can extend up to 20" when pulled down. Maximum weight supportable in the retracted position is 15 ounces.

Hose Retractor 15-61-6006

More Allied Innovations

Suction Regulators



Timeter® Flowmeters



Gomco® Suction Equipment



B&F Oxygen Therapy Disposables



LSP Bag Mask Resuscitators



Brass Regulators & Oxygen Cylinders



Allied Healthcare Products, Inc. • 1720 Sublette Avenue • St. Louis, Missouri 63110
314/771-2400 • 800/444-3954 • Fax: 314/771-0650 • www.alliedhpi.com
International offices in Europe, Asia and South America

U.S.

Toll Free: 800/444-3940
Fax: 800/477-7701

Canada & Mexico

Toll Free: 800/446-0552
Fax: 800/246-6201

International

Telephone: 314/268-1683
Fax: 314/771-5183