



Columbia Dentoform® Teaching Solutions

DentalEZ® Equipment

NevinLabs™ Workstations

RAMVAC® Utility

StarDental® Instruments

Affordable, LubeFree Dry Vacuums



- LubeFree and waterless
- Small footprint
- Open frame design for ease of maintenance
- Stackable units for easy expansion
- Auto shut down procedure to ensure longevity of rotary vanes
- 5 Year Warranty

RAMVAC® Utility Solutions



RAMVAC products are manufactured in an ISO 13485:2003 Certified facility.

SOF_7.2.3-01-201



Benefits for Replacing Water Ring Systems:

- Dry system reduces water and sewage costs
- No traps to clean
- Consistent vacuum levels when one or more HVE's are opened
- Easy installation
- Able to adapt easily to most existing plumbing (3/4" – 2")

Benefits for New Dry Vacuum Applications:

- Separation tank and E1 control are all integrated into complete unit
- Can use flexible PVC
- Light weight
- Economically priced

Badger Specs (Part # BADG1.2.60.1)	
Voltage	208-230 VAC, Single Phase
Full Load	6.7 A
Frequency	60 Hz
Horsepower	1.5 HP
Maximum Vacuum	15 in Hg
Preset Vacuum Level	7 in Hg
Dimension	18.7" W x 25.1" D x 22.7" H
Tank Capacity	5 gallons
Tank Material	Polyethylene
Unit Weight	152 lbs
Shipping Weight	195 lbs

Includes: LubeFree rotary vane vacuum pump, 5 gallon integrated separation tank, and E1 control.

Meets the following requirements:

- UL 60601-1 (2nd Edition)
- CSA C22.2 No. 601.1
- NFPA 99c Level 3 Vacuum System

