


## MED TOUCH HMI ALARM PANELS







- Master / Area Alarms
- Combination Alarms
- Up to 128 signal points in a single panel
- Bright, easy to read LCD displays
- Text/e-mail notification of alarm events
- Enter text instructions for staff under alarm condition.
- Wireless models available
- Best Warranty: 5 Year parts / 2 Year Labor



- Easy setup & install
- Transducers may be installed remotely or in back box
- User programmable Hi/Lo Setpoints
- Easily input custom text instructions





Feature	Medical Gas Manufacturer			
	<u> </u>	Beacon Medaes	Amico	Ohio Medical
Touch Screen	<b>1</b>	<b>1</b>	X	<b>*</b>
Area Alarm	<b>/</b>	1	<b>/</b>	<b>*</b>
Maximum Number of Gas/Area Boards	(16)	8	8	8
Master Alarm	<b>*</b>	<b>*</b>	<b>/</b>	X
Maximum Number of Master Signal Points	(128)	60 (one panel)	30	X
Combination Alarm	<b>1</b>	<b>*</b>	X	X
Transducer Master	1	火	X	X
4 – 20 mA Master	<b>/</b>	<b>-</b>	X	X
Emergency Instructions	1	<b>-</b>	X	X
Text Notifications	1	1	X	X
Ethernet Option	1	<b>1</b>	<b>/</b>	<b>-</b>
Website	1	1	<b>1</b>	K
Event Log	1	1	X	X
Modbus	1	X	X	<b>-</b>
BACnet	1	1	Separate box – expensive	X
Last Event History	1	X	X	K
Repeater on Master Alarm Signals	1	X	X	X
Wireless Master LCD Alarm	1	X	X	X
Integra Touch LCD Area Alarm (Combination Area Alarm & Zone Valve Box in One)	✓	×	×	X

For more information, contact your local Powerex distributor or Powerex sales representative today!



150 Production Dr • Harrison, OH 45030 powerexinc.com • (T)888.769.7979 • (F)513.367.3125

